

Medtronic Bedside Remote Monitoring

- Please keep your bedside monitor always plugged in to an area that receives adequate Wi-Fi and cellular signal.
- How to send a manual transmission (only by the request of your physician):
 - Sit within 6 feet of your monitor, and make sure it is securely plugged in on both ends of the wire.
 - Press and hold the center button for 3 seconds and release. A **GREEN** light will spin around the button.
 - Wait near the monitor until all the bars below the button fill to solid **GREEN**.
 - When the light around the button is fully solid **GREEN**, no further action is needed.
 - If you see any **RED** color on the monitor, call Medtronic support during business hours for assistance.
- Segment of light below the center button and the three-bar icon is illuminated **YELLOW**: Your monitor cannot connect to a network. Try to send a manual transmission later.
- **GREEN** segment of light below the center button blinks three times: Your implanted device cannot connect to the monitor. Try to send a manual transmission later.
- Segment of light below the center button and a person with phone illuminated **YELLOW**: Your monitor failed a 'self-check' and may need to be replaced. Call Medtronic support.

Medtronic Support: 866-470-7709

